



Digital Microscope

ShuttlePix

P-400Rv



A new era of microscopy: "Shuttle style."

ShuttlePix

ShuttlePix, a revolutionary portable digital microscope

The ShuttlePix digital microscope is an all-new, one-of-a-kind digital microscope that can be used as a docked device, for high magnification and various measurements, or as a portable digital camera, for capturing any sample sizes.



Zoom Camera Head

- O Cordless setup enabling image capture in the lab or in the field (Battery-powered with built-in illumination and SD card slot)
- ODigital camera-like ease of use

ShuttlePix

• 4-segment LED ring light

High-resolution image capture

o 20x optical zoom

• Scene mode

Motorized Stand


- Extended Depth-of-Focus (EDF) image capture
- Intuitive stylus and icon-driven operation using dedicated Touch Panel Monitor (advanced image capturing, simple measurement, etc.)
- Dedicated PC software expands operation and possibilities (image processing, 3D/color height maps, etc.)

Simple Stand

- All-in-one design not needing a power source (simple reflection stand)
- Provides diascopic illumination imaging capabilities (simple diascopic LED stand)
- Ouick view feature ensures fast image uploads to your PC
- Easy navigation with the Remote Control NEW


Experience the high performance by Nikon's high quality optics and advanced technologies


ShuttlePix STYLE



$\equiv 20x$ optical zoom


When used on the motorized focusing stand and touch panel monitor, this magnification increases to between 20x and 400x. ShuttlePix provides ample magnification range that ensures you never have to switch or purchase different lenses for different magnifications. Magnification information is also linked to ShuttlePix's scale and simple measurement functions.








Ensure optimal settings for image capture through five types of scene mode: wafer/IC chip, metal, printed circuit board, flat panel display, and high-reflection sample.



Printed circuit board



High-reflection sample NEW (polyethylene form)

High NA / High definition / Wide field of view **Maximum optical performance**

Nikon's proprietary optics achieve precise observation and imaging with NA up to 0.2 (at 400x magnification) and 20mmdiagonal wide field of view (at 20x magnification). Changing of Resolution Preferred Mode and Depth-of-Focus Preferred Mode is also possible.










Newly-developed design for bright, even illumination

4-segment LED ring light

ShuttlePix's new illumination technology achieves consistent brightness at all levels of magnification. Capture shaded images as well, through split-half illumination switchable among top, bottom, left, and right.





Molded plastic

From shooting large samples to use in the field Capable of shooting any sample size

The Handy can observe samples of all sizes, such as inspections of art displays to automotives. Use of the stand allows for abrasion observation of electronic and aircraft parts. With the built-in battery, cable-less operation is possible.
















3 "Shuttle Styles" To Freely Choose From

High-resolution, all-focused image capturing is possible without PC

with MOTORIZED STAND



■ Motorized Focusing Stand plus Touch Panel Monitor


ShuttlePix is equipped with a vertical-movement Motorized Focusing Stand and a 17", 1280×1024 color LCD Touch Panel Monitor. Through the intuitive operation of touching icons or using the screen stylus, precise image capturing and simple measurement are now possible.





EDF images can be captured simply by operating the button located on the stand.








Quick 3D view NEW

Captured images are simultaneously transmitted to the Shuttle Pix Editor, where 3D images can be viewed.





Perform measurements without a PC

Simple measurements

The touch panel monitor makes it possible to perform a truly wide range of simple measurements.



Measurement function

Distance between 2 points, point-to-line distance, distance between center points of 2 circles, angle,

Annotation

Count marking, text input, pen drawing, straight line, airs, grid, XY scale, XY


Stage options for every application

Stage lineup

Image capturing for large samples (up to 75mm×50mm×148mm) is possible with the Zoom Camera Head and Motorized Focusing Stand. Select from three dedicated Stages matched to the observation sample.









Capture digital microscopy images at any location

ZOOM CAMERA HEAD

Image capture with the ease of a digital camera
Zoom Camera Head

Image capturing is simple, with the use of the capture switch (trigger). The lightweight and ergonomic camera head allows for easy handling for all users. Voice recording is now available as a new function.

■ Scale display

(cross hairs, H, grid)

The scale display feature allows for

confirmation of a sample or defect's





Achieving optimal image capture
Best shot selection mode

When capturing images, up to 10 continuous frames are recorded and only the optimal shot is kept. This guarantees capturing sharp images even at high magnifications.





(scheduled release in August 2013)

During Preview Mode, images will be automatically captured at the optimum focus point.



High magnification image capturing and simple measurements anywhere anytime

with SIMPLE STAND

■ Stands tailored to the purpose of use **■** Simple reflection/ Diascopic LED stand

An all-in-one compact reflector stand with a battery operated Zoom Camera Head and a Diascopic LED Stand capable of diascopic observation are both available.





Automatically upload your images to PC Quick view

The free, dedicated software ShuttlePix Editor (see p.6) can automatically upload images taken with the ShuttlePix to your PC. This enables users to immediately scroll through thumbnail images and perform simple measurements on selected images right after image capturing.



Camera Head.

Wide range of features using the dedicated software

with PC / SOFTWARE

Free download and registration of software **■** ShuttlePix Editor



Conveniently outputs simple measurements and 3D EDF images or cross-section displays directly into Excel via dedicated ShuttlePix Editor* software. Software download from the Nikon website and user registration are free.

*Compatible with Windows XP and Windows 7 (Quick view feature only compatible with Windows® 7)

Speed capture makes it easy to measure right after shooting an image **Interface**

Zoom camera head set: Use the zoom head camera's SD card to upload saved images to your PC.


Motorized focusing stand set: In addition to SD card, use the USB memory connected to the stand or USB cable to upload saved images to your PC.


Plain stand set: In addition to the SD card and USB memory, use the Quick view feature to upload images to your PC automatically.



3D display of EDF images and height data 3D display

Use a 3D bird's-eye view to display EDF images and height data taken with ShuttlePix. Rotation, zoom in/out, scale display, colorbased heightmaps, and other image display operations are available.







3D display

Color heightmap

Display and simple measurement of height data cross-section Cross-section display and simple measurement

Display cross sections at specified positions based on height data embedded in EDF images. Perform simple measurements of the cross section including height, angle, and width, with measurement data displayable and recordable in tabular form.





Original image

Cross section graph

Various measuring performance **Simple measurement**

Add comments and markers to key measurements. Measurement results can be output in tabular form.













Angle Rule Horizontal ... 0.4299


Vertical Ru... 0.6639






0.4299







See http://www.nikon.com/products/instruments/lineup/digital-microscope/shuttlepix/index.htm


for details and download information regarding ShuttlePix Editor software. *Download of ShuttlePix Editor is limited to registered users of the product.




Zoom Camera Head









Zoom Camera Head (P-400Rv)	
Effective pixels	Approx. 2 megapixels
CCD	CCD: 1/1.8" color CCD, total pixels approx. 2.11 megapixels
	Frame rate: 28fps (selectable 28fps (800×600) or 15fps (1600×1200) when
	connected to Motorized Focusing Stand)
Optics	Magnification: 3× to 60× (magnification on built-in 2.7" monitor),
	20× to 400× (magnification on dedicated 17" monitor), optical zoom ratio = 20:1
Working distance	
FOV	Maximum diagonal field of view 20mm (16mm x 12mm)
Illumination	Light source: White LED
	Illumination method: Episcopic illumination from around the objective lens
	Illumination area: ø20mm, 4-segment ring LED (top/bottom/left/right)
Recording	Storage media: SD memory card, SDHC memory card (max. 16GB)
	(selectable USB memory or FTP when connected to Motorized Focusing Stand)
	Image file format: TIFF (non-compressed), JPEG (3 compression levels)
	Audio file format: WAV (voice memo)
	Recording pixels: 2M (1600×1200), 0.5M (800×600)
Shooting mode	Scene mode (standard, wafer/ IC chip, metal ceramic, circuit board, FPD,
	high-reflection sample*), BSS (Best Shot Select), timer (2 seconds fixed), interval,
_	4 custom settings, Focus assist capture function* (scheduled release in August 2013)
Exposure	Photometry method: Average photometry/peak-hold photometry
	Exposure control: Program AE/shutter preferred/manual exposure
	Exposure compensation: -2EV to +2EV in 1/3EV steps,
	camera gain and shutter speed can be set (manual exposure)
	AE lock function
Aperture	Resolution preferred mode/Depth-of-focus preferred mode

Image	Image quality adjustment: Saturation/hue/contrast/sharpness/color effect
compensation	Shading correction: Factory setting (switchable On/Off can be switched)
White balance	Manual setting (adjustable red/blue gain can be adjusted)
LCD Monitor	2.7" TFT color LCD, turned off automatically when connected to Motorized Focusing Stand
Image playback	Full-frame view, thumbnail view (9 frames), zoom view (scrollable)
Image deletion	Quick delete, select image delete, folder delete, card format
Video output	NTSC/PAL
Connector	Video output, capture input (ø3.5mm stereo mini jack), dedicated connector for stand, external remote controller input
Language	Japanese/English
Power supply	Li-ion Rechargeable Battery/AC adapter/Motorized Focusing Stand (when connected to Motorized Focusing Stand)
	Battery: P-Li01 P-Li-ion Rechargeable Battery
	AC adapter: AC adapter EA1050E-120 (optional)
	Auto power save: 30sec./1min./5min./10min./20min.
Power consumption	24VA
Battery operating time	Approx. 90 minutes (battery life at maximum power consumption with maximum LED brightness)
Charging time	Approx. 4 hours (when no charge remains)
Tripod socket hole	1/4-inch (ISO 1222)
Dimensions	Approx. 96(W)×228(D)×142(H)mm
Weight	Approx. 1kg (excluding battery and SD card)
Usage environment	Temperature: 0 to +40°C; Humidity: 60% RH max. (no condensation)

^{*}Zoom camera head, "P-400R", can also be used by downloading the firmware from the Nikon Website

Motorized Focusing Stand (P-MFSC)

	• , ,
Stroke	Z axis stroke: 150mm (upward 148mm, downward 2mm), upper and lower limit adjustable
Stage	3×2 Stage/Sliding Stage/Tilting Stage
Image Edit	EDF: Still image display, 3D image display (with ShuttlePix Editor)
	Halation prevention: Live display (max. 5fps, 1280×960, 800×600)
	High dynamic range: Still image display (1280×960/800×600)
	Calibration: Zoom magnification conversion/user-registerd calibration
	Measurement function: Distance between 2 points, point-to-line distance, distance between center points of 2 circles, angle, circle, area, pitch
	Annotation: Count marking, text input, pen drawing, straight line, scale indication, cross-hairs, grid, XY scale, XY measurement
Direct Printing	Supported printer: PictBridge printer

Touch Panel Monitor (P-TPM)

Display size	17.0" (display area 333.9mm×270.3mm)
Resolution	SXGA (1280×1024)
Connector	Image input: Digital input: DVI-D SXGA, video input: NTSC composite
	USB host: USB2.0 A connector×3
	USB device: USB2.0 B connector×1

Plain Stand C (P-PSC)

Dimensions	Approx. 180(W)×265(D)×325.5(H)mm
Stroke	Z axis stroke: 100mm (upward 98mm, downward 2mm),
Size of stage glass	ø90mm
Weight	Approx. 2.5kg

Connector	Camera head IF: for connection to P-400Rv
	Video output: DVI-I
	LAN: 10Base-T/100Base-TX
	USB host: USB2.0 A connector×2 (for connection to USB mouse/ USB memory/ dedicated monitor)
	USB device: USB2.0 B connector×1 (for connection to PC)
Language	Japanese/English
Power supply	Built-in, AC100-240V 50/60Hz
	Auto power save: OFF/1 to 99min.
Power consumption	140VA
Dimensions	Approx. 280(W)×404(D)×450(H)mm
Weight	Approx. 9kg
Usage environment	Temperature: 0 to +40°C; Humidity: 60% RH max. (no condensation)

Power supply	AC100-240V 50/60Hz
Power consumption	140VA
Dimensions	Approx. 390(W)×180(D)×371(H)mm
Weight	Approx. 7kg
Usage environment	Temperature: 0 to +40°C; Humidity: 60% RH max. (no condensation)

LED Stand (P-LEDS)

Dimensions	Approx. 200(W)×325(D)×350.5(H)mm
Stroke	Z axis stroke: 100mm (upward 98mm, downward 2mm),
Size of stage glass	ø90mm
Weight	Approx. 4.5kg

http://www.nikon.com/products/instruments/lineup/digital-microscope/shuttlepix/index.htm

See the URL above for product details and download information regarding ShuttlePix Editor dedicated software. *Download of ShuttlePix Editor is limited to registered users of the product.



Received the iF Design Award's Product Design Award 2012

Received the Good Design Award 2011

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. July 2013 ©2010-2013 NIKON CORPORATION

N.B. Export of the products* in this catarog is controlled under the Japanese Foreign Exchange and Foreign Trade Law. Appropriate export procedures shall be required in case of export from Japan.
*Products: Hardware and its technical information (including software)



TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING THE EQUIPMENT.



NIKON CORPORATION

Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan phone: +81-3-6433-3701 fax: +81-3-6433-3784 http://www.nikon.com/instruments/





NIKON METROLOGY, INC.

12701 Grand River Avenue, Brighton, MI 48116 U.S.A. phone: +1-810-220-4360 fax: +1-810-220-4300 E-mail: Sales.US.NM@nikon.com

http://www.nikonmetrology.com/

NIKON METROLOGY EUROPE NV

Geldenaaksebaan 329, 3001 Leuven, Belgium phone: +32-16-74-01-00 fax: +32-16-74-01-03 Email: Sales.Europe.NM@nikon.com http://www.nikonmetrology.com/

NIKON INSTRUMENTS (SHANGHAI) CO., LTD.

CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060 (Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026 (Guangzhou branch) phone: +86-20-3882-0550 fax: +86-20-3882-0580

NIKON INSTRUMENTS KOREA CO., LTD.

KOREA phone: +82-2-2186-8400 fax: +82-2-555-4415

NIKON SINGAPORE PTE LTD.

SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668

NIKON MALAYSIA SDN. BHD.

MALAYSIA phone: +60-3-7809-3688 fax: +60-3-7809-3633

PT. NIKON INDONESIA

INDONESIA phone: +62-21-574-6262 fax: +62-21-574-6363 Nikon Sales (Thailand) Co., Ltd.

THAILAND phone: +66-2633-5100 fax: +66-2633-5191

NIKON INDIA PRIVATE LIMITED

INDIA phone: +91-124-4688500 fax: +91-124-4688527

NIKON CANADA INC.

CANADA phone: +1-905-602-9676 fax: +1-905-602-9953

NIKON INSTRUMENTS S.p.A. ITALY phone: +39-055-300-96-01 fax: +39-055-30-09-93

NIKON METROLOGY UK LTD.

UNITED KINGDOM phone: +44-1332-811-349 fax: +44-1332-639-881 E-mail: Sales.UK.NM@nikon.com

NIKON METROLOGY SARL

FRANCE phone: +33-1-60-86-09-76 fax: +33-1-60-86-57-35

E-mail: Sales.France.NM@nikon.com

NIKON METROLOGY GMBH

GERMANY phone: +49-6023-91733-0 fax: +49-6023-91733-229 E-mail: Sales.Germany.NM@nikon.com