



Nikon Digital Sight Series New Lineup

A new system for imaging: the DS-Fi3, a high resolution and sensitivity general purpose color camera has been added to the Nikon Digital Sight series. The DS-Fi3 can be connected to a PC, or the new compact tablet-style DS-L4.



Microscope Camera DS-Fi3 Nikon High-resolutio 5.9

Color







Tubular adenoma, HE staining (Objective: CFI Plan Apochromat λ 4x) Photos courtesy of:Dr. Yasunori Ohta, Department of Pathology, IMSUT Hospital, Institute of Medical Science, The University of Tokyo

DS-FI3



Breast cancer, FISH method (Objective: CFI Plan Apochromat 100x Oil) Photos courtesy of: Hironao Kusakari, Diagnostic Pathology, St. Marianna University Hospital

High sensitivity, low noise

Quantum efficiency and read noise have been greatly improved, providing better capability for acquisition of fluorescent images with better signal-to-noise ratios than before.



eft image: Uterine cervix Pap. Staining

(Objective: CFI Plan Apochromat λ60x Oil)

Pathology, Yokohama City University Hospital

(Objective: CFI Plan Achromat NCG 40x) Photos courtesy of: Clinical Laboratory

rtment, Yokohama City University Hospital

Right image: Bone marrow

Photos courtesy of: Kazuhiro Mita, Department of

High-speed live display

Fast USB3.0 data transfer means fast, smooth live updating of images for finding samples or focusing, even at full resolution.



Superior color reproduction

Nikon is well-known for outstanding and lifelike color reproduction, and developing superior algorithms for creating results that look like the actual samples. These algorithms are used in all of the color cameras in the digital sight lineup.





Camera Control

The DS-FI3 interfaces with PC computers via a USB3.0 interface directly to the camera head, and uses NIS-Elements series software for image acquisition.



Compact, easy-to-use tablet-type microscope camera control unit.

DS-Fi3 can be optionally connected to the DS-L4 tablet-style control unit, eliminating the need and space requirements of a desktop PC. DS-L4 has a large number of built-in functions for measurement and annotations, and has built-in security for network connectivity.

Tablet-type camera control unit

Large, 10.1 inch, touch-screen 1920 × 1200 pixel display: The DS-Fi3 can be set and operated simply and easily through the tablet by touch, or by connecting Bluetooth accessories such as a keyboard or mouse.

User Interface for naturally simple operation

The camera control menu uses recognizable and intuitive icons.. Frequently used icons are in two rows, and the display space for live images and photographed images is large and prominently displayed.



Live menu/Replay menu

Scene mode

When connected to biological, industrial, or stereoscopic microscopes equipped with motorized hardware units and observation mode sensors, it is possible to both control the microscope and detect its observation mode state. Storing the objective lens information is convenient when making measurements.

Biological Scene Mo	de	
Brightfield	HE	ELISA
Industrial Scene Mod	le	
Wafer/IC	Z Meta	I, Ceramic/Plastic
Circuit board	🕎 Flat F	Panel Display

Integration with microscopes

Optimal imaging parameters for the microscope's light source, (LED or halogen), each sample type, and observation method can easily be set through the icons. A choice of three modes for biological

and four modes for industrial imaging are available, and up to seven custom modes with freely configurable shooting parameters can be set.

	8	162		-		
200 200 200	5. 1				1	1
24 01K 0.000		10		Į R		
21 Let 0.000		125		18		
÷		0				

Used with ECLIPSE Ni-E

A wide variety of tools

The DS-L4 enables easy measurements directly on images, with input of lines and comments. These can also be written and saved with the image, and measurement data can also be output.



Measurement (2 point distance)

Tighten security

McAfee embedded control with White List method is preinstalled for the virus measurement.

The program which is not registered at White List cannot be launched so that the virus cannot be activated.

Only registered users are able to use by implementing user registration.

Those security protect the important images.

Two Large Sensor high resolution 16.25-megapixel CMOS image sensors for microscopy

Two Nikon FX-format CMOS image sensor cameras join the Digital Sight series of microscope digital cameras: the DS-Ri2 color digital camera and the DS-Qi2 monochrome digital camera.

High pixel density and large field of view coupled with USB3.0 high speed data transfer offer fast frame rates and high resolution images with these CMOS image sensors.



Florida State University

DS-Ri2

16.25 megapixel (not interpolated) and accurate color rendition are features that make the DS-Ri2 an excellent choice for recreating color images as they eyes see them.



Large Format **CMOS** image sensors

Nikon manufactures CMOS image sensors and imaging technologies for professional DSLR cameras, and now has optimized our sensors for microscopy



Mikan

DS-Ria

DS-Qi2

High pixel density, high sensitivity and low noise are key features of the DS-Qi2 monochrome camera.



Pig kidney epithelial cells expressing GFP-EB3 tubulin Sample courtesy of: Michael Davidson, National High Magnetic Field Laboratory,

The tissues of the liver, HE staining (Objective: CFI Plan Apochromat λ 10x) Photos courtesy of: Kazuhiro Muraoka, Photography Division, Tokyo Women's Medical University

Fast, one-shot capture of ultra-high resolution color images.

Microscope Camera





Mouse cerebellum sagittal section, HE staining (Objective: CFI Plan Apochromat λ 4x)

High-resolution images

16.25-megapixel CMOS image sensors for astonishing image quality

The DS series enables one-shot instantaneous capture and fast storage of images with resolution as high as 4908 x 3264 pixels, without pixel shifting or pixel stepping.

This pixel density is ideally suited for photomicrography of ultra-fine structures or patterns in biological or industrial samples, at low or high magnifications.

Photography with the natural colors seen through the microscope

Nikon is a leader in development of algorithms for reproducing color just as the eyes see it

The DS models' image processing engine is based on extensive data accumulated over many years of developing microscope color digital cameras, resulting in perfect reproduction of the colors your eyes see in the microscope.





Pancreatic cancer cell, NGFR immunostaining*1 (Objective: CFI Plan Apochromat λ 40x)

Human glomerulus of kidney, Azan stain*2 (Objective: CFI Plan Apochromat λ 40x)

*1,*2 Photos courtesy of: Dr. Atsushi Furuhata and Noriyoshi Sueyoshi, Assistant General Manager, Laboratory of morphology and image analysis, BioMedical Research Center, Juntendo University Graduate School of Medicine





Sensitivity settings that span the range from ISO200 to ISO12800 allow the capture of vivid fluorescent color images.

> Transgenic C. elegans expressing venus in the head neurons and EGFP in the body wall muscles. Photos courtesy of: Drs. Keiko Gengyo-Ando and Junichi Nakai, Saitama University Brain Science Institute

Capture Low light fluorescence and Large Fields of View





Indian Muntjac Deer Skin Fibroblast Cells, Cytoskeletal F-actin labeled with Alexa Fluor 488 Sample courtesy of: Michael Davidson and Florida State University

High sensitivity

Detects even faint fluorescent signals

7.3µm pixels, high quantum efficiency, and very low read noise allow the DS-Qi2 to read in even faint fluorescent signals.



Low noise

Acquires dim fluorescent signals with ultra-low noise

Both 2.2 electrons read noise coupled with a large full-well capacity and 0.6 electrons dark current allow the acquisition of 14bit fluorescence images with very little noise.

> LLC-PK1 cells expressing GFP-EB3 tubulin with low noise. Large linear full well capacity allows acquiring both the brightest and dimmest areas in a single capture. Sample courtesy of: Michael Davidson, National High Magnetic Field Laboratory, Florida State University

Excellent linearity

Reliable quantitative analysis made possible

With a linearity error of ±1%, the DS-Qi2 is a superb tool for measuring intensities in fluorescence samples, including timebased intensity measurement and ratiometric measurement.

High frame rate

Fast focusing, even with fluorescent images

With a high-sensitivity CMOS image sensor and USB 3.0-based data transfer, the DS-Qi2 enables high-speed live imaging and image capture at up to 45 fps (1636×1088 pixels).







Fluorescent time-lapse imaging through integration with NIS-Elements software

imaging applications.



Rat primary culture neuron Dendron labeled with MAP-2(Red) and Actin(cytoskeleton) labeled with Phalloidin (Green)

LLC-PK1 cells expressing GFP-EB3 tubulin (green) and H2B-labeled histones (red) illustrating the large field of view of the DS-Qi2 camera. Sample courtesy of: Michael Davidson, National High Magnetic Field Laboratory, Florida State University





Integration with the comprehensive imaging software series

Nikon uses the NIS-Elements series as control software. NIS-Elements allows functions from basic imaging to control of the microscope and peripheral devices to be performed, as well as the measurement, analysis, and management of acquired images. Four basic packages and a variety of optional modules are available to suit every application and objective. * See the NIS-Elements Catalog for details.





The bundled free package offers functions for the display of scale on live images, full-screen display, and more. The simple operation screen makes shooting easy.



with measurement and report creation

functions. It enables general microscopic

image acquisition in fields from biomedical to

industrial, and is expandable through optional

added features such as EDF and databases.



The research package enables the construction of advanced image acquisition systems, including multidimensional imaging (up to 4 dimensions for Br, 6 dimensions for Ar), through integration with systemized microscopes. Sets equipped with a rich range of image processing and analysis functions are available for every application.

Compatible OS: Windows® 7 Pro 32/64bit

* Nikon provides confirmed compatible PCs with up-to-date specifications. Contact Nikon for details.

Multichannel (multi color) Ar Br

NIS-Elements can acquire full bit depth multi-color images, combining multiple fluorescence wavelengths and different illumination methods (DIC, phase contrast etc.), while offering independently scalable channels.







All-color merged image

Z-series Ar Br D

Through motorized focus control, NIS-Elements reconstructs and renders 3D images from multiple Z-axis planes.







Multi-dimensional Image Display Ar Br

NIS-Elements displays time lapse, multi-channel, multiple X, Y, Z positions in an intuitive layout, which allows for automatic playback and the ability to select subsections of the data to be saved as a new file





22. call*/ 12.ants 546 22. Works Bril

HDR (High Dynamic Range) image acquisition



HDR creates an image with appropriate brightness in both the dark and bright regions in a sample by combining multiple images acquired with different exposure settings. It is also possible to create HDR image using multiple captured images.



Captures both areas **1** and **2** with optimal exposure

EDF (Extended Depth of Focus)

Option Ar Br D

Creates a single, all-in-focus image from images of differing focus. Such images can now be created by simply turning the focus knob.





Selects the in-focus area and produces one all-in-focus image

Image stitching (Large Image)

Option Ar Br D

Stitches together images from multiple fields of view during shooting to create an image with wide field of view. Images already acquired can also be stitched together.



Manual measurement and image annotation Ar Br D



Manual Measurement allows easy measurement of length and area by drawing lines or an object directly on the image. The results can be attached to the image, and also exported as text or to an Excel spreadsheet.



Auto measurement (Object Counting) Ar Br Option D

Performs binarization on images using previously set thresholds to measure the number, area, brightness, etc. of identified objects.





Grain size analysis

Detects and measures grains in one and two phase samples according to JIS G0551 or ASTM E112-96/E1382-97 standards.





Cast iron analysis

Option Ar Br D

Option Ar Br D

Detects, measures and classifies graphite content as well as ferrite content in graphite-corrected samples according to JIS G5502 or ASTM A247-06 standards



Dimensions

DS-L4



0.0

E

265.5

10.2



DS-Ri2/DS-Qi2

DS-L4 cradle AC adapter



System Diagram



Specifications

Model name	DS-Fi3	DS-Ri2	DS-Qi2	
Image sensor	1/1.8 inch Color CMOS image sensor Size: 6.91 × 4.92 mm	Nikon FX-format Color CMOS image sensor Size: 36.0 × 23.9 mm	Nikon FX-format Monochrome CMOS image sensor Size: 36.0 × 23.9 mm	
Recordable pixels	All pixels: 2880 × 2048 2 Vertical and 2 horizontal pixels average: 1440 × 1024	All pixels: 4908 × 3264 3 × 3 pixels average: 1636 × 1088		
Lens mount	C-mount	F-mount		
Cooling method		_	Electronic cooling	
ISO sensitivity (recommended exposure index)	Standard: equivalent to ISO 50 (Selectable from ISO 50 to ISO 3200 equivalent)	Standard: equivalent to ISO 200 (Selectable from ISO 200 to ISO 12800 equivalent)	Standard: equivalent to ISO 800 (Selectable from ISO 800 to ISO 51200 equivalent)	
Quantum efficiency		_	77%	
Full well Capacity		_		
Readout noise		_	2.2e (- typ.)	
Dark current		_		
Live display mode*1 (maximum fps)	All pixels (2880 × 2048): 15 fps 2 Vertical and 2 horizontal pixels average (1440 × 1024): 30 fps	All pixels (4908 × 3264): 6 fps 3 × 3 pixels average (1636 × 1088): 45 fps		
Exposure time	100 µsec ~ 30 sec	100 µsec ~ 120 sec		
Photometry mode	Average photometry: Average intensity within	n the photometry area Peak photometry:	Maximum intensity within the photometry area	
Exposure control	One-time automatic exposure: Exposure tim Continuous automatic exposure: Automatic e Manual exposure: Exposure time and gain s	e is adjusted automatically for one-time w exposure adjustment is performed continu ettings are made manually	ithin the optimum range for the camera Jously to keep the exposure within the camera	
Exposure correction	±1EV Step:1/6EV		Average metering: -1 EV ~ +1/2 EV Peak hold metering: -1 EV ~ ±0 EV	
Interface	USB3.0 (connect with PC, DS-L4*2) × 1, External trigger × 1		USB3.0 (connect with PC) \times 1, External trigger \times 1	
Power supply	AC100-240V 50Hz/60Hz			
Power consumption	4.8 W	13W	24W	
Dimensions	100(W) × 66(D) × 65(H)mm	105(W) × 134(D) × 153(H)mm		
Weight	400g (approx.)	1200g (approx.)		
Operating environment	0-40°C, 60% RH max. (without condensation)		0-30 °C, 80% RH max. 30-40°C, 60% RH max. (without condensation	

Microscope Camera Control Unit -----

Model name	
Connectable cameras*2	DS-Fi3
Live image	FULL (resolution emphasized), FAST (display frame rate
Exposure control	Program AE/ Focus AE/Manual AE With AE lock functio
Brightness adjustment	Exposure mode AE: Exposure compensation adjustmen
Exposure metering	Average metering, Peak hold metering
Exposure metering area	Position/size variable
White balance	One-push operation
Image correction	Tone, sharpness, black level, hue, chroma, R/B adjustme
Scene mode	LED/halogen, Biological/industrial/asbestos/standard
Recording format	Colorspace: sRGB Still image: Tiff/Jpeg/DICOM, movie
Saving destination	Internal drive (flash memory 64GB, approx. 28GB availa
Measurement/drawing/ scale	Measurement target: Point-to-point distance, perpendicu Measurementunit calibration registration: Auto registratio Manual registration calibration: Manual calibration/Optic Drawing: Text, line, arrow, pen, marker, scale bar Scale: Cross, glid line, X scale (cross scale), XY scale
Microscope control*4	Biological microscope: Ni-E/Ni-U/Ci-E Industrial micros Stereo microscope: SMZ25/SMZ18/SMZ1270i
Supported language	English, Japanese
Security	Anti-virus: McAfee Embedded Control is preinstalled. Pr User login method: With DS-L4 user registration, login is
LCD display	10.1-inch wide TFT LCD display (1920×1200 WUXGA)
Interface	USB 3.0 host port ×1, USB 2.0 port ×4, DisplayPort, mic Bluetooth, microSD card slot
Power supply	AC100-240V 50Hz/60Hz
Power consumption	65W
Dimensions	Main body: $265 \times 188 \times 10$ mm, Including extended cra
Weight	Main body: 630 g (approx.), Including extended cradle:
Operating environment	5-35°C, 20-80% RH max. (without condensation)

*1: Maximum frame rate depends on exposure time.

*3: Shading correction will be avilable by updating a software version in January 2017.

DS-L4

emphasized) Format: RGB 24bits

nt, Exposure mode manual AE: Exposure time or gain adjustment

nent, shading correction*3, special effect Custom: Up to 7 types can be registered e: AVI

able) /microSD/USB memory/SMB file server

ular line length, angle, circle, distance, between center points of two circles, area ion calibration by objective information setting (seven types registerable) cal calibration by entering objective magnification (14types)

oscope: LV150NA/LV100DA-U/LV100NDA

rograms that are started by the white list method are controlled. possible by entering user ID and password

croHDMI, LAN (IEEE 802.3 10/100/1000Base-TX, IEEE 802.11 a/b/g/n),

adle: 265 × 201 × 107 mm 1060 g (approx.)

*2: DS-Ri2 will be able to connect to DS-L4 in January 2017. *4: Please contact Nikon for details.

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. September 2016 ©2004-2016 NIKON CORPORATION

N.B. Export of the products⁺ in this catalog is controlled under the Japanese Foreign Exchange and Foreign Trade Law. Appropriate export procedures shall be required in case of export from Japan. *Products: Hardware and its technical information (including software)

WARNING TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING THE EQUIPMENT.



Shinagawa Intercity Tower C, 2-15-3, Konan, Minato-ku, Tokyo 108-6290, Japan Microscope Solutions Business Unit, phone: +81-3-6433-3705 fax: +81-3-6433-3785 Industrial Metrology Business Unit, phone: +81-3-6433-3701 fax: +81-3-6433-3784 http://www.nikon.com/instruments/



Nikon

NIKON INSTRUMENTS INC. 1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A. phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only) fax: +1-631-547-0306

http://www.nikoninstruments.com/

NIKON METROLOGY, INC. 12701 Grand River Avenue, Brighton, MI 48116 U.S.A. phone: +1-810-220-4360 fax: +1-810-220-4300

E-mail: Sales.US.NM@nikon.com http://www.nikonmetrology.com/ NIKON INSTRUMENTS EUROPE B.V.

Tripolis 100, Burgerweeshuispad 101, 1076 ER Amsterdam, The Netherlands phone: +31-20-7099-000 fax: +31-20-7099-298 http://www.nikoninstruments.eu/

NIKON METROLOGY EUROPE NV

Geldenaaksebaan 329, 3001 Leuven, Belgium phone: +32-16-74-01-00 fax: +32-16-74-01-03 Email: Sales.Europe.NM@nikon.com http://www.nikonmetrology.com/

NIKON INSTRUMENTS (SHANGHAI) CO., LTD.

CHINA phone: +86-21-6841-2050 fax: +86-21-6841-2060 (Beijing branch) phone: +86-10-5831-2028 fax: +86-10-5831-2026 (Guangzhou branch) phone: +86-20-3882-0550 fax: +86-20-3882-0580 **NIKON INSTRUMENTS KOREA CO., LTD.** KOREA phone: +82-2-2186-8400 fax: +82-2-555-4415

NIKON SINGAPORE PTE LTD.

SINGAPORE phone: +65-6559-3651 fax: +65-6559-3668 NIKON MALAYSIA SDN. BHD. MALAYSIA phone: +60-3-7809-3688 fax: +60-3-7809-3633 PT. NIKON INDONESIA INDONESIA phone: +62-21-574-6262 fax: +62-21-574-6363 Nikon Sales (Thailand) Co., Ltd. THAILAND phone: +66-2633-5100 fax: +66-2633-5191 NIKON INDIA PRIVATE LIMITED INDIA phone: +91-124-4688500 fax: +91-124-4688527 NIKON CANADA INC. CANADA phone: +1-905-602-9676 fax: +1-905-602-9953 NIKON INSTRUMENTS S.p.A. ITALY phone: +39-55-300-96-01 fax: +39-55-30-09-93 NIKON AG SWITZERLAND phone: +41-43-277-28-67 fax: +41-43-277-28-61 NIKON GMBH AUSTRIA AUSTRIA phone: +43-1-972-6111-00 fax: +43-1-972-6111-40 NIKON BELUX BELGIUM phone: +32-2-705-56-65 fax: +32-2-726-66-45 NIKON UK LTD.

UNITED KINGDOM phone: +44-208-247-1717 fax: +44-208-541-4584 NIKON METROLOGY UK LTD.

UNITED KINGDOM phone: +44-1332-811-349 fax: +44-1332-639-881 E-mail: Sales.UK.NM@nikon.com





MIX

Paper from

FSC[•] C013657

IS0 14001 Certified for NIKON CORPORATION

NIKON FRANCE S.A.S.

FRANCE phone: +33-1-4516-45-16 fax: +33-1-4516-45-55 **NIKON METROLOGY SARL** FRANCE phone: +33-1-60-86-09-76 fax: +33-1-60-86-57-35

E-mail: Sales.France.NM@nikon.com NIKON GMBH

GERMANY phone: +49-211-941-42-20 fax:+49-211-941-43-22 NIKON METROLOGY GMBH GERMANY phone: +49-6023-91733-0 fax: +49-6023-91733-229 E-mail: Sales.Germany.NM®nikon.com