

EXTEC Abrasive Powders are formulated for use in distilled water mixtures and available in standard particle sizes and grades. Aluminum Oxide (Al₂O₃) and Silicon Carbide (SiC) Powders are recommended for fine grinding and rough polishing.

EXTEC Magnesium Oxide (MgO) is recommended when fine polishing wrought aluminum alloys and magnesium based metals and alloys.

EXTEC Iron Oxide (Fe₂O₃) is recommended to polish soft ores, glass, lapidary, and soft metals.

EXTEC Chromic Oxide (Cr₂O₃) is recommended when Aluminum Oxide contamination is not permissible.



| CATALOG NUMBER | DESCRIPTION | UNIT |
|--|-------------|---------------|
| EXTEC Aluminum Oxide (Al₂O₃) Powder | | |
| 17012 | 3.0 µm | 1 lb (.45 kg) |
| 17016 | 3.0 µm | 5 lb (2.3 kg) |
| 17020 | 5.0 µm | 1 lb (.45 kg) |
| 17024 | 5.0 µm | 5 lb (2.3 kg) |
| 17028 | 9.5 µm | 1 lb (.45 kg) |
| 17032 | 9.5 µm | 5 lb (2.3 kg) |
| 17036 | 12.5 µm | 1 lb (.45 kg) |
| 17040 | 12.5 µm | 5 lb (2.3 kg) |
| 17044 | 14.5 µm | 1 lb (.45 kg) |
| 17048 | 14.5 µm | 5 lb (2.3 kg) |
| 17052 | 17.5 µm | 1 lb (.45 kg) |
| 17056 | 17.5 µm | 5 lb (2.3 kg) |
| 17060 | 20.0 µm | 1 lb (.45 kg) |
| 17064 | 20.0 µm | 5 lb (2.3 kg) |
| 17076 | 25.0 µm | 1 lb (.45 kg) |
| 17080 | 25.0 µm | 5 lb (2.3 kg) |
| 16836 | 240 Grit | 1 lb (.45 kg) |
| 16838 | 240 Grit | 5 lb (2.3 kg) |
| 16846 | 400 Grit | 1 lb (.45 kg) |
| 16848 | 400 Grit | 5 lb (2.3 kg) |
| 16852 | 600 Grit | 1 lb (.45 kg) |
| 16854 | 600 Grit | 5 lb (2.3 kg) |

| CATALOG NUMBER | DESCRIPTION | UNIT |
|--|---|---------------|
| EXTEC Silicon Carbide (SiC) Abrasive Powder | | |
| 16802 | 120 Grit | 1 lb (.45 kg) |
| 16804 | 120 Grit | 5 lb (2.3 kg) |
| 16806 | 240 Grit | 1 lb (.45 kg) |
| 16808 | 240 Grit | 5 lb (2.3 kg) |
| 16812 | 320 Grit | 1 lb (.45 kg) |
| 16814 | 320 Grit | 5 lb (2.3 kg) |
| 16816 | 400 Grit | 1 lb (.45 kg) |
| 16818 | 400 Grit | 5 lb (2.3 kg) |
| 16822 | 600 Grit | 1 lb (.45 kg) |
| 16824 | 600 Grit | 5 lb (2.3 kg) |
| 16826 | 1000 Grit | 1 lb (.45 kg) |
| 16828 | 1000 Grit | 5 lb (2.3 kg) |
| EXTEC Magnesium Oxide (MgO) Abrasive Powder | | |
| 16910 | EXTEC Magnesium Oxide (MgO) Abrasive Powder | 1 lb (.45 kg) |
| EXTEC Levigated Aluminum Oxide (Al₂O₃) Abrasive Powder 5 µm | | |
| 16903 | EXTEC Levigated Aluminum Oxide (Al ₂ O ₃) Abrasive Powder 5 µm | 1 lb (.45 kg) |
| 16905 | | 5 lb (2.3 kg) |
| EXTEC Iron Oxide (Fe₂O₃) Abrasive Powder 3 µm | | |
| 16915 | EXTEC Iron Oxide (Fe ₂ O ₃) Abrasive Powder 3 µm | 1 lb (.45 kg) |
| 16918 | | 5 lb (2.3 kg) |
| EXTEC Chromic Oxide (Cr₂O₃) Abrasive Powder 1-5 µm | | |
| 16943 | EXTEC Chromic Oxide (Cr ₂ O ₃) Abrasive Powder 1-5 µm | 1 lb (.45 kg) |
| 16945 | | 5 lb (2.3 kg) |