

EXTEC Coated Abrasives are specifically designed for wet or dry (W/D) material sample preparation and feature a durable, heavy weight or poly backing. Available in belt, disc, strip, roll and sheet form; Pressure Sensitive Adhesive (PSA) backed in disc or strip form.

EXTEC Silicon Carbide (SiC) Abrasives are ideal for rough and fine grinding and polishing.

EXTEC Aluminum Oxide (Al<sub>2</sub>O<sub>3</sub>) Abrasives are ideal for rough and fine grinding.

EXTEC Zirconia Aluminum Abrasives are heavy duty and aggressive.



CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Belt SiC W/D 4" x 36" (102 mm x 914 mm)</b>		
13745	50 Grit	Pkg/10
13750	80 Grit	Pkg/10
13755	120 Grit	Pkg/10
13760	180 Grit	Pkg/10
13765	240 Grit	Pkg/10
13770	320 Grit	Pkg/10
13775	400 Grit	Pkg/10
13780	600 Grit	Pkg/10
<b>EXTEC Abrasive Belt Al<sub>2</sub>O<sub>3</sub> W/D 4" x 36" (102 mm x 914 mm)</b>		
13705	50 Grit	Pkg/10
13710	80 Grit	Pkg/10
13715	120 Grit	Pkg/10
13720	180 Grit	Pkg/10
13725	240 Grit	Pkg/10
<b>EXTEC Abrasive Belt Zirconia Alumina W/D 4" x 36" (102 mm x 914 mm)</b>		
13785	60 Grit	Pkg/5
13790	120 Grit	Pkg/5
13795	180 Grit	Pkg/5
13800	220 Grit	Pkg/5

CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Disc SiC W/D PSA 2" (51 mm)</b>		
15251	60 Grit	Pkg/100
15252	120 Grit	Pkg/100
15253	180 Grit	Pkg/100
15254	240 Grit	Pkg/100
15255	320 Grit	Pkg/100
15256	400 Grit	Pkg/100
15257	600 Grit	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D PSA 2 7/8" (73 mm)</b>		
15262	120 Grit	Pkg/100
15263	180 Grit	Pkg/100
15264	240 Grit	Pkg/100
15266	320 Grit	Pkg/100
15267	400 Grit	Pkg/100
15268	600 Grit	Pkg/100
15270	10 Each 180,240 320,400,600 Grit	Pkg/50 Assortment

CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Disc SiC W/D 8" (203 mm)</b>		
15001	60 Grit	Pkg/100
15002	120 Grit	Pkg/100
15003	180 Grit	Pkg/100
15004	240 Grit	Pkg/100
15005	320 Grit	Pkg/100
15006	400 Grit	Pkg/100
15007	600 Grit	Pkg/100
15008	1200 Grit	Pkg/100
15010	800/2400 Grit	Pkg/100
15012	1200/2400 Grit	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D Poly Back 8" (203 mm)</b>		
15015	5 µm	Pkg/100
15017	1 µm	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D PSA 8" (203 mm)</b>		
15032	120 Grit	Pkg/100
15033	180 Grit	Pkg/100
15034	240 Grit	Pkg/100
15035	320 Grit	Pkg/100
15036	400 Grit	Pkg/100
15037	600 Grit	Pkg/100
15038	1200 Grit	Pkg/100
15040	800/2400 Grit	Pkg/100
15042	1200/4000 Grit	Pkg/100
15043	10 Each 240, 320, 400, 600 Grit	Pkg/40 Assortment
<b>EXTEC Abrasive Disc SiC W/D 8" with 1/2" Center Hole (203 mm x 12.7 mm)</b>		
15171	60 Grit	Pkg/100
15172	120 Grit	Pkg/100
15173	180 Grit	Pkg/100
15174	240 Grit	Pkg/100
15176	320 Grit	Pkg/100
15177	400 Grit	Pkg/100
15178	600 Grit	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D PSA 8" with 2 1/2" Center Hole (203 mm x 64 mm)</b>		
15222	120 Grit	Pkg/100
15223	180 Grit	Pkg/100
15224	240 Grit	Pkg/100
15226	320 Grit	Pkg/100
15227	400 Grit	Pkg/100
15228	600 Grit	Pkg/100
<b>EXTEC Abrasive Disc Zirconia Alumina W/D PSA 8" (203 mm)</b>		
15044	60 Grit	Pkg/25
15045	120 Grit	Pkg/25
15046	180 Grit	Pkg/25
15047	220 Grit	Pkg/25

CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Disc SiC W/D 12" (305 mm)</b>		
15081	60 Grit	Pkg/100
15082	120 Grit	Pkg/100
15083	180 Grit	Pkg/100
15084	240 Grit	Pkg/100
15085	320 Grit	Pkg/100
15086	400 Grit	Pkg/100
15087	600 Grit	Pkg/100
15088	1200 Grit	Pkg/100
15090	800/2400 Grit	Pkg/100
15091	1200/4000 Grit	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D Poly Back 12" (305 mm)</b>		
15092	5 µm	Pkg/100
15093	1 µm	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D 12" with 1/2" Center Hole (305 mm x 12.7 mm)</b>		
15094	60 Grit	Pkg/100
15095	120 Grit	Pkg/100
15096	180 Grit	Pkg/100
15097	240 Grit	Pkg/100
15098	320 Grit	Pkg/100
15099	400 Grit	Pkg/100
15100	600 Grit	Pkg/100
<b>EXTEC Abrasive Disc SiC W/D PSA 12" (305 mm)</b>		
15141	60 Grit	Pkg/100
15142	120 Grit	Pkg/100
15143	180 Grit	Pkg/100
15144	240 Grit	Pkg/100
15145	320 Grit	Pkg/100
15146	400 Grit	Pkg/100
15147	600 Grit	Pkg/100
15148	1200 Grit	Pkg/100
15150	800/2400 Grit	Pkg/100
15152	1200/4000 Grit	Pkg/100
<b>EXTEC Abrasive Disc Zirconia Alumina W/D PSA 12" (305 mm)</b>		
15153	60 Grit	Pkg/25
15154	120 Grit	Pkg/25
15155	180 Grit	Pkg/25
15156	220 Grit	Pkg/25



CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Strip SiC W/D PSA 3" x 11" (76 mm x 279 mm)</b>		
15244	240 Grit	Pkg/100
15246	320 Grit	Pkg/100
15247	400 Grit	Pkg/100
15248	600 Grit	Pkg/100
15249	25 Each 240, 320, 400, 600 Grit	Pkg/100
<b>EXTEC Abrasive Sheet SiC W/D 9" x 11" (229 mm x 279 mm)</b>		
15181	60 Grit	Pkg/100
15182	120 Grit	Pkg/100
15183	180 Grit	Pkg/100
15184	240 Grit	Pkg/100
15186	320 Grit	Pkg/100
15187	400 Grit	Pkg/100
15188	600 Grit	Pkg/100

CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Disc SiC W/D PSA 16" (406 mm)</b>		
15380	60 Grit	Pkg/100
15381	120 Grit	Pkg/100
15382	180 Grit	Pkg/100
15383	240 Grit	Pkg/100
15385	320 Grit	Pkg/100
15386	400 Grit	Pkg/100
15387	600 Grit	Pkg/100
<b>EXTEC Abrasive Roll SiC W/D 3" x 75' (76 mm x 23 m)</b>		
15190	240 Grit	Each
15191	320 Grit	Each
15192	400 Grit	Each
15193	600 Grit	Each
<b>EXTEC Abrasive Roll SiC W/D 3-7/16" x 75' (88 mm x 23 m)</b>		
15194	240 Grit	Each
15196	320 Grit	Each
15197	400 Grit	Each
15198	600 Grit	Each

CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Abrasive Disc Al<sub>2</sub>O<sub>3</sub> W/D Poly Back PSA 4" (102 mm)</b>		
20985	12 μm	Pkg/50
20986	9 μm	Pkg/50
20987	3 μm	Pkg/50
20988	1 μm	Pkg/50
20989	0.3 μm	Pkg/50
<b>EXTEC Abrasive Disc Al<sub>2</sub>O<sub>3</sub> W/D Poly Back PSA 8" (203 mm)</b>		
21045	12 μm	Pkg/50
21046	9 μm	Pkg/50
21047	3 μm	Pkg/50
21048	1 μm	Pkg/50
21049	0.3 μm	Pkg/50
<b>EXTEC Abrasive Disc Al<sub>2</sub>O<sub>3</sub> W/D Poly Back PSA 12" (305 mm)</b>		
21005	12 μm	Pkg/50
21006	9 μm	Pkg/50
21007	3 μm	Pkg/50
21008	1 μm	Pkg/50
21009	0.3 μm	Pkg/50
<b>EXTEC Abrasive Disc Al<sub>2</sub>O<sub>3</sub> W/D Poly Back PSA 9" x 11" (229 mm x 279 mm)</b>		
21025	12 μm	Pkg/25
21026	9 μm	Pkg/25
21027	3 μm	Pkg/25
21028	1 μm	Pkg/25
21029	0.3 μm	Pkg/25