

Thermosetting EXTEC Edge Retention Epoxy Molding Compounds exhibit excellent flow, minimum shrinkage, chemical resistance, high hardness and superior adhesion.

EXTEC Ultra Fine Edge Retention Epoxy Molding Compound is recommended for encapsulating particularly intricate specimens.

EXTEC Diallyl Phthalate Mineral Powder provides moderate shrinkage and chemical resistance. Diallyl Phthalate Glass is formulated with short glass for a harder surface.

Conductive Thermosetting EXTEC Mounting Compounds are available in carbon and copper compounds.

Transparent/clear and fragile applications require EXTEC Transparent Thermoplastic Mounting Powder.

EXTEC Thermosetting Phenolic Compression Powders are available in Black, Red, Green, and Mottled Green/Brown. 1 1/4" and 1 1/2" Black, Red, and Green Preforms are easy to handle and accelerate molding time.



CATALOG NUMBER	DESCRIPTION	UNIT
14327	EXTEC Hydraulic Fluid	qt (.95 l)
14330	EXTEC Mounting Gloves	Pair
14335 ◇	EXTEC Silicon Mold Release Spray	Can 6 oz (.17 kg)
14340	EXTEC Metal Thermometer 0-250° C.	Each
14505	EXTEC Edge Retention Epoxy Molding Compound (Black)	5 lb (2.3 kg)
14500	EXTEC Ultra Fine Edge Retention Epoxy Compound (Black)	5 lb (2.3 kg)
14497	EXTEC Diallyl Phthalate Short Glass Mounting Powder (Blue)	5 lb (2.3 kg)
14494	EXTEC Blue Diallyl Phthalate Mineral Mounting Powder (Blue)	5 lb (2.3 kg)
14509	EXTEC Carbon Conductive Mounting Compound (Black)	1 lb (.45 kg)
14503	EXTEC Copper Conductive Mounting Compound	5 lb (2.3 kg)
14515	EXTEC Transparent Thermoplastic Mounting Powder	5 lb (2.3 kg)

◇ Special Packaging Required

CATALOG NUMBER	DESCRIPTION	UNIT
<b>EXTEC Black Phenolic Mounting Powder</b>		
14350	Powder	5 lb (2.3 kg)
14352	Powder	25 lb (11.5 kg)
<b>EXTEC 1 1/4" (32 mm) Black Phenolic Preform</b>		
14370	Preform	100/pc
14372	Preform	500/pc
14375	Preform	1000/pc
<b>EXTEC 1 1/2" (38 mm) Black Phenolic Preform</b>		
14380	Preform	100/pc
14382	Preform	500/pc
14385	Preform	1000/pc
<b>EXTEC Green Phenolic Powder</b>		
14450	Powder	5 lb (2.3 kg)
14452	Powder	25 lb (11.5 kg)
<b>EXTEC 1 1/4" (32 mm) Green Phenolic Preform</b>		
14470	Preform	100/pc
14472	Preform	500/pc
14475	Preform	1000/pc
<b>EXTEC 1 1/2" (38 mm) Green Phenolic Preform</b>		
14480	Preform	100/pc
14482	Preform	500/pc
14485	Preform	1000/pc
<b>EXTEC Red Phenolic Powder</b>		
14400	Powder	5 lb (2.3 kg)
14402	Powder	25 lb (11.5 kg)
<b>EXTEC 1 1/4" (32 mm) Red Phenolic Preform</b>		
14420	Preform	100/pc
14422	Preform	500/pc
14425	Preform	1000/pc
<b>EXTEC 1 1/2" (38 mm) Red Phenolic Preform</b>		
14430	Preform	100/pc
14432	Preform	500/pc
14435	Preform	1000/pc
14486	EXTEC Green/Brown Mottled Phenolic Mounting Powder	5 lb (2.3 kg)