


EXTEC exceptional polishing cloths and pads come in all standard sizes and textures for every laboratory application. Plain and adhesive backs. Special-size orders accepted.



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Alphagam™ Cloth		
A medium nap cloth with synthetic Rayon fibers bonded to a woven cotton cloth. For final polishing of ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, and pcb's. Recommended for use with EXTEC Diamond Products and EXTEC Polishing Powders and Suspensions.		
17238	For 8" (203 mm) Wheel	Pkg/10
17290	For 250 mm 10" Wheel	Pkg/10
17292	For 12" (305 mm) Wheel	Pkg/10
17294	For 14" (356 mm) Wheel	Pkg/10
EXTEC Alphagam™ Cloth, PSA		
Same Versatile cloth as above with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant.		
17239	For 2" (50 mm) Wheel PSA	Pkg/20
17240	For 2 7/8" (73 mm) Wheel PSA	Pkg/20
17254	For 4" (100 mm) Wheel PSA	Pkg/20
17256	For 6" (152 mm) Wheel PSA	Pkg/10
17258	For 8" (203 mm) Wheel PSA	Pkg/10
17259	For 8" (203 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17260	For 250 mm 10" Wheel PSA	Pkg/10
17262	For 12" (305 mm) Wheel PSA	Pkg/10
17282	For 12" (305 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17284	For 14" (356 mm) Wheel PSA	Pkg/10
17286	For 16" (406 mm) Wheel PSA	Pkg/10




CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Ultragam™ Cloth PSA		
A fine fiber, high density, short nap cloth that results in maximum flatness with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For final polishing of ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, and pcb's. Recommended for use with EXTEC Diamond Products and EXTEC Polishing Powders and Suspensions.		
17287	For 8" (203 mm) Wheel PSA	Pkg/10
17288	For 250 mm 10" Wheel PSA	Pkg/10
17289	For 12" (305 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Multi-Cloth™ PSA A short nap, synthetic velvet material that features a thin durable back with Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant and provides minimal relief and edge rounding. For final polishing of ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, and pcb's. Recommended for use with EXTEC Diamond Products and EXTEC Polishing Powders and Suspensions.		
17602	For 2" (50 mm) Wheel PSA	Pkg/20
17607	For 2 7/8" (73 mm) Wheel PSA	Pkg/20
17606	For 8" (203 mm) Wheel PSA	Pkg/10
17610	For 8" (203 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17608	For 250 mm 10" Wheel PSA	Pkg/10
17611	For 12" (305 mm) Wheel PSA	Pkg/10
17612	For 12" (305 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17614	For 14" (356 mm) Wheel PSA	Pkg/10




CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Plano™ Pad PSA A napless chemotextile material that is oil and water resistant with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For final polishing of ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, pcb's, minerals and glass. Recommended for use with EXTEC Diamond Products.		
17500	For 2" (50 mm) Wheel PSA	Pkg/20
17502	For 2 7/8" (73 mm) Wheel PSA	Pkg/20
17504	For 4" (100 mm) Wheel PSA	Pkg/20
17506	For 6" (152 mm) Wheel PSA	Pkg/10
17518	For 8" (203 mm) Wheel PSA	Pkg/10
17519	For 8" (203 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17520	For 250 mm 10" Wheel PSA	Pkg/10
17522	For 12" (305 mm) Wheel PSA	Pkg/10
17542	For 12" (305 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17544	For 14" (356 mm) Wheel PSA	Pkg/10
17545	For 16" (406 mm) Wheel PSA	Pkg/10
EXTEC Plano Perforated™ Pad PSA Similar to EXTEC Plano with Perforations with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For final polishing of ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, pcb's, minerals and glass. Recommended for use with EXTEC Diamond Products.		
17557	For 6" (152 mm) Wheel PSA	Pkg/20
17558	For 8" (203 mm) Wheel PSA	Pkg/10
17559	For 8" (203 mm) Wheel with 2 1/2 (6.4 cm) Center Hole PSA	Pkg/10
17560	For 250 mm 10" Wheel PSA	Pkg/10
17562	For 12" (305 mm) Wheel PSA	Pkg/10
17582	For 12" (305 mm) Wheel with 2 1/2 (6.4 cm) Center Hole PSA	Pkg/10
17584	For 14" (356 mm) Wheel PSA	Pkg/10
17586	For 16" (406 mm) Wheel PSA	Pkg/10




CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Duraplan Plano™ Pad PSA A LONG LASTING napless Chemotextile material that is oil and water resistant. Features Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For intermediate polishing of hard ferrous metals, ceramics, carbides and thermal spray coatings. Recommended for use with EXTEC Diamond Products.		
17588	For 2 7/8" (73 mm) Wheel PSA	Pkg/20
17590	For 6" (152 mm) Wheel PSA	Pkg/10
17592	For 8" (203 mm) Wheel PSA	Pkg/10
17594	For 250 mm 10" Wheel PSA	Pkg/10
17596	For 12" (305 mm) Wheel PSA	Pkg/10
17598	For 14" (356 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Duraplan Ultra™ Cloth PSA A napless, hard-woven durable synthetic cloth that is extremely aggressive with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For intermediate polishing of hard ferrous metals, ceramics, carbides and thermal spray coatings. Recommended for use with EXTEC Diamond Products.		
17192	For 8" (203 mm) Wheel PSA	Pkg/10
17194	For 250 mm 10" Wheel PSA	Pkg/10
17196	For 12" (305 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Optigam™ Cloth PSA A closely woven, napless synthetic cloth bonded to a poly barrier with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For general purpose final polishing of hard ferrous metals, ceramics, carbides, and thermal spray coatings. Recommended for use with EXTEC Diamond Products.		
17300	For 6" (152 mm) Wheel PSA	Pkg/10
17302	For 8" (203 mm) Wheel PSA	Pkg/10
17304	For 250 mm 10" Wheel PSA	Pkg/10
17306	For 12" (305 mm) Wheel PSA	Pkg/10
17308	For 14" (356 mm) Wheel PSA	Pkg/10
17310	For 16" (406 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC S-Plan™ Cloth PSA A napless, exceptionally tight woven synthetic cloth featuring an extremely smooth surface, free from roughness with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For final polishing ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, and pcb's. Recommended for use with EXTEC Diamond Products.		
17400	For 8" (203 mm) Wheel PSA	Pkg/10
17402	For 250 mm 10" Wheel PSA	Pkg/10
17404	For 12" (305 mm) Wheel PSA	Pkg/10




CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Chemic-Cloth™ PSA A napless, chemically resistant synthetic cloth for maximum flatness with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant for chemomechanical polishing techniques. For final polishing of hard metals, soft metals, cast iron, composites, thermal spray coatings, electronics, pcb's, ceramics, minerals, glass and polymers. Recommended for use with EXTEC Diamond Products and EXTEC Polishing Powders and Suspensions.		
17668	For 8" (203 mm) Wheel PSA	Pkg/10
17670	For 8" (203 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17671	For 250 mm 10" Wheel PSA	Pkg/10
17672	For 12" (305 mm) Wheel PSA	Pkg/10
17674	For 14" (356 mm) Wheel PSA	Pkg/10
17676	For 16" (406 mm) Wheel PSA	Pkg/10




CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Poly Cloth™ PSA A short nap, resilient flocked fiber material coated to a polyester film backing with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For final polishing of ferrous and nonferrous metals, composites, polymers, cast irons, ceramics, carbides, thermal spray coatings, electronics, and pcb's. Recommended for us with EXTEC Diamond Products and EXTEC Polishing Powders and Suspensions.		
17618	For 8" (203 mm) Wheel PSA	Pkg/10
17624	For 250 mm 10" Wheel PSA	Pkg/10
17625	For 12" (305 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Durawool™ Cloth PSA A high quality durable woven wool cloth with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For final polishing of ferrous and nonferrous metals, composites, thermal spray coatings, polymers, minerals and glass. Recommended for use with EXTEC Diamond Products and EXTEC Polishing Powders and Suspensions.		
17599	For 8" (203 mm) Wheel PSA	Pkg/5
17600	For 250 mm 10" Wheel PSA	Pkg/5
17601	For 12" (305 mm) Wheel PSA	Pkg/5




CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Red Felt Cloth PSA A buffed, 100% fine virgin wool cloth with uniform thickness designed for long life with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For general rough polishing of steel and other ferrous metals. Recommended for use with EXTEC Polishing Powders and Suspensions.		
17348	For 8" (203 mm) Wheel PSA	Pkg/10
17350	For 250 mm 10" Wheel PSA	Pkg/10
17352	For 12" (305 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Cot-Cloth™ A napless, woven cotton cloth recommended for rough polishing or softer materials. Recommended for use with EXTEC Diamond Products, Extec Polishing Powders and Suspensions.		
17208	For 8" (203 mm) Wheel	Pkg/10
17210	For 250 mm 10" Wheel	Pkg/10
17212	For 12" (305 mm) Wheel	Pkg/10
17216	For 14" (356 mm) Wheel	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Synthetic Velvet Cloth A high pile, premium synthetic velvet with a woven base for final polishing of soft materials. Recommended for the final polishing of soft materials.		
17418	For 8" (203 mm) Wheel	Pkg/10
17420	For 250 mm 10" Wheel	Pkg/10
17422	For 12" (305 mm) Wheel	Pkg/10
17424	For 14" (356 mm) Wheel	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Billiard Cloth PSA A sheared, 100% virgin wool, high quality cloth with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For general rough polishing. Recommended for use with EXTEC Polishing Powders and Suspensions.		
17325	For 2" (50 mm) Wheel PSA	Pkg/20
17330	For 8" (203 mm) Wheel PSA	Pkg/10
17332	For 250 mm 10" Wheel PSA	Pkg/10
17334	For 12" (305 mm) Wheel PSA	Pkg/10
17323	For 14" (356 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Nylon Cloth PSA A napless, long lasting cloth with a Pressure Sensitive Adhesive (PSA) backing that is oil, water & alcohol resistant. For rough and final polishing. Recommended for use with EXTEC Polishing Powders and Suspensions.		
17156	For 2 7/8" (73 mm) Wheel PSA	Pkg/20
17166	For 6" (152 mm) Wheel PSA	Pkg/10
17168	For 8" (203 mm) Wheel PSA	Pkg/10
17178	For 8" (203 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17170	For 250 mm 10" Wheel PSA	Pkg/10
17172	For 12" (305 mm) Wheel PSA	Pkg/10
17182	For 12" (305 mm) Wheel PSA with 2 1/2 (6.4 cm) Center Hole	Pkg/10
17185	For 14" (356 mm) Wheel PSA	Pkg/10
17190	For 16" (406 mm) Wheel PSA	Pkg/10



CATALOG NUMBER	DESCRIPTION	UNIT
EXTEC Silk Cloth A 100% natural silk cloth recommended for final polishing. Recommended for use with EXTEC Diamond Products.		
17378	For 8" (203 mm) Wheel	Pkg/10
17380	For 250 mm 10" Wheel	Pkg/10
17382	For 12" (305 mm) Wheel	Pkg/10
17384	For 14" (356 mm) Wheel	Pkg/10

